

## Aspire COVE Pressure Relief Chair



- The Aspire Cover chair helps to improve the postural, pressure care and comfort outcomes for users that require medium to high level support.
- Users can be effortlessly moved between rooms at home, or in care environments. This offers improved inclusivity and social interaction.
- **Anatomical Legrest Movement:** As the legrest raises, it extends outwards tracking knee joint range of motion and maximises lower limb support. When combined with correct seat depth adjustment, taller users maintain access under tables promoting social inclusiveness and dining independence.
- **Infinite Tilt Adjustment:** The gas assisted tilting mechanism provides a smooth and intuitive motion to maximise user comfort and confidence. By tilting up to 40 degrees, sitting pressure is redistributed from vulnerable bony prominences to the lower back and assists in maintaining postural alignment.
- **Configurable Pressure Support Cushioning:** The integrated seat cushion combines layered memory foam with a segmented hybrid air cell that offers convenient pressure adjustment to carers and maintenance staff. Reinflation is as simple as opening and closing a valve.
- **Ease of Movement:** High performance TENTE Levina castors have been chosen to reduce push/pull forces during transportation. This improves workplace health and safety outcomes for care staff and family members.
- Directional locking castors at the front of the chair assist during transit, particularly on uneven flooring or carpets with pronounced grains.
- **Infection Control:** The sophisticated medical grade upholstery has been chosen to perform under modern infection control protocols. Removable seating elements and smooth wipeable surfaces simplify cleaning processes between users. The biostatic patient surface prevents bacteria and virus ingress.
- **Market Leading Upholstery:** The patient surface is genuine DARTEX made in the UK. This breathable highly elastic surface maximises skin microclimate and pressure redistribution.

- External surfaces are medical grade coated fabrics from Warwick Australia.
- Hardy yet stylish patterns offer an appealing aesthetic and blend seamlessly in with modern Aged Care interior designs.
- The layered seat structure accommodates after-market pressure care products such as ROHO cushions, without compromising chair performance or infection control outcomes (purchased separately).
- Transitional pressure management armrests.
- Ergonomic carer push handles.
- Tool free adjustments.
- Aspire Cove Chair Accessories available here. {insert link to CHP137}

## Warranty Information

- Fabric: 1 year
- Frame: 3 Years

| Code         | Display Name                                                     | Total Length | Total Width | Total Height | Safe Working Load | Colour              | Seat Depth | Seat Width | Seat To Floor Height | Seat To Footrest | Backrest Height |
|--------------|------------------------------------------------------------------|--------------|-------------|--------------|-------------------|---------------------|------------|------------|----------------------|------------------|-----------------|
| UK-CHP198960 | Aspire COVE Pressure Relief Chair - Medium - Grey/Andreas Pumice | 33" - 64"    | 31"         | 4'2"         | 25 st             | Grey/Andreas Pumice | 20" - 23"  | 22"        | 21"                  | 21"              | 28" - 30"       |
| UK-CHP198965 | Aspire COVE Pressure Relief Chair - Medium - Navy/Bronx Steel    | 33" - 64"    | 31"         | 4'2"         | 25 st             | Navy/Bronx Steel    | 20" - 23"  | 22"        | 21"                  | 21"              | 28" - 30"       |
| UK-CHP198970 | Aspire COVE Pressure Relief Chair - Large - Grey/Andreas Pumice  | 33" - 66"    | 32"         | 4'3"         | 25 st             | Grey/Andreas Pumice | 20" - 23"  | 24"        | 21"                  | 21"              | 28" - 30"       |
| UK-CHP198975 | Aspire COVE Pressure Relief Chair                                | 33" - 66"    | 32"         | 4'3"         | 25 st             | Navy/Bronx          | 20" - 23"  | 24"        | 21"                  | 21"              | 28" - 30"       |

| Code | Display Name               | Total Length | Total Width | Total Height | Safe Working Load | Colour | Seat Depth | Seat Width | Seat To Floor Height | Seat To Footrest | Backrest Height |
|------|----------------------------|--------------|-------------|--------------|-------------------|--------|------------|------------|----------------------|------------------|-----------------|
|      | - Large - Navy/Bronx Steel |              |             |              |                   | Steel  |            |            |                      |                  |                 |